





Chirana Progress
Vrbovská cesta 17
921 01 Piešťany
Slovakia

tel.: +421 33 7956211
e-mail: chirana@chirana-progress.sk
www.chirana-progress.sk



180



Useful capacity	280 l
Water Pump	433 l/min / 0.36 MPa
Dimensions	2500 x 960 x 920-942 mm
Power input	1800 VA
Power supply	3 x 400 V / 50 (60) Hz

Underwater **massage**



24

180



Useful capacity	280 l
Water Pump	433 l/min / 0.36 MPa
Dimensions	2500 x 960 x 920-942 mm
Power input	3500 VA
Power supply	3 x 400 V / 50 (60) Hz

Underwater and whirling **massage**



40

180



3

Useful capacity	280 l
Water Pump	433 l/min / 0.48 MPa
Dimensions	2500 x 960 x 920-942 mm
Power input	3500 VA
Power supply	3 x 400 V / 50 (60) Hz

Underwater and whirling **massage**



80/100

180



2/4

Useful capacity	280 l
Water Pump	433 l/min / 0.48 MPa
Dimensions	2500 x 960 x 920-942 mm
Power input	3500 VA
Power supply	3 x 400 V / 50 (60) Hz

Automatic underwater and whirling **massage**



174

180



5

Useful capacity	260 l
Water Pump	433 l/min / 0.48 MPa
Dimensions	2500 x 960 x 920-942 mm
Power input	3500 VA
Power supply	3 x 400 V / 50 (60) Hz

Automatic underwater and whirling **massage**

aquaviva®

Ocean® Economy

Ocean® Standard

Ocean® Forte

Ocean® Forte PC

Ocean® de Luxe



Models **Laguna Compact**
A more compact version of the Laguna. Shorter by 25 centimeters and with flat bottom.

Available in the same versions as regular sized Laguna bath tubs (Plus, Bubble, Plus Bubble), except for multizone versions.

Models Laguna Compact	
Useful capacity	200 l
Dimensions	1950 x 870 x 790-815 mm
Models Laguna	
Useful capacity	200 l
Dimensions	2200 x 870 x 790-815 mm

Balneological **bath tub**



25

2

Useful capacity	200 l
Water Pump	433 l/min / 0.36 MPa
Dimensions	2200 x 870 x 790-815 mm
Power input	3500 VA
Power supply	3 x 400 V / 50 (60) Hz

Combined balneological and whirling **bath tub**



180

Useful capacity	200 l
Water Pump	433 l/min / 0.36 MPa
Dimensions	2200 x 870 x 790-815 mm
Power input	450 VA
Power supply	230 V / 50 (60) Hz

Combined balneological and whirling **bath tub**



25

180

2

Useful capacity	200 l
Water Pump	433 l/min / 0.36 MPa
Dimensions	2200 x 870 x 790-815 mm
Power input	3950 VA
Power supply	3 x 400 V / 50 (60) Hz

Combined balneological and whirling **bath tub**



124

180

5

Useful capacity	200 l
Water Pump	450 l/min / 0.35 MPa
Dimensions	2200 x 870 x 790-815 mm
Power input	3500 VA
Power supply	3 x 400 V / 50 (60) Hz

Automatic underwater and whirling **massage**



Water jets



Air jets



174

Number of jets



Underwater massage



Standard equipment



Optional equipment

5

Max no. of massage zones



Electrotherapy

Optional equipment, depending on the model, includes:

- RGB chromato therapy
- Air massage
- Automatic filling
- Pulsation
- Sine / wave / ymphatic massage
- Electromagnetic incrust stopper
- Ozone therapy
- Magneto therapy
- Connection to and from a recycling system
- RGB backlit bath tub
- RGB backlit LCD
- Automatic disinfection
- Log file - data about bath tub operation

Note:

- Some models are available in 230 V version. Contact us for more details.

All photos shown are for illustration purpose only. Parameters of the product can change depending on the optional equipment.

More detailed information can be found at:
www.chirana-progress.sk

Laguna® Compact

Laguna® Plus

Laguna® Bubble

Laguna® Plus Bubble

Laguna® Tornado

Explanatory notes



Useful capacity 20 l
Dimensions 2500 x 960 x 920 mm
Power input 50 VA
Power supply 230 V / 50 (60) Hz

Whole body galvanic / diadynamic **bath**

Electra® CG



Useful capacity 700 l
Water Pump 233 l/min / 0,48 MPa
Dimensions 2700 x 1930 x 980 mm
Power input 2400 VA
Power supply 3 x 400 V / 50 (60) Hz

Underwater and whirling **massage**

Hubbard bath



Useful capacity 370 l
Water Pump 133 l/min / 0,57 MPa
Dimensions 2600 x 1000 x 1200 mm
Power input 1800 VA
Power supply 3 x 400 V / 50 (60) Hz

Underwater **massage**

VOD 56®



Useful capacity 300 l
Water Pump 133 l/min / 0,57 MPa
Dimensions 1935 x 1500 x 760 mm
Power input 30 VA
Power supply 230 V / 50 (60) Hz

Balneological / birthing **bath tub**

Hubbard Plus



Water supply pressure 0,3 - 0,5 MPa
Dimensions 500 x 1000 x 1120 mm
Distance to patient 4 - 5 m

Spraying **desk**

Niagara®, Niagara® Plus



Capacity 14 - 21 l
Operating temp. in peat pack holders 40 - 45 °C
Operating temperature for paraffin 56 - 58 °C
Temperature control range 0 - 90 °C
Power supply 230 V / 50 Hz
Power input 1500 - 2 x 1500 VA

Paraffin and peat pack **warmer**

Tep®



Current types: a) Galvanic DC current
b) Sine
c) Faraday
d) Triangle
e) Saw

Additional settings for “b” to “e” currents:
MF monophase,
DF diphasic,
CP short period,
LP long period

Useful capacity arms / legs 13 / 24,5 l
Dimensions 1280 x 950 x 975 mm
Power input 30 VA
Power supply 230 V / 50 (60) Hz

4-chamber galvanic / diadynamic **bath**

Electra®




Kneipp: Alternating baths of the limbs in empty tanks with contrasting hot and cold water showers..

Kneipp: Bath in water with steadily increasing temperature.

Useful capacity arms / legs 13 / 24,5 l
Dimensions 1280 x 950 x 975 mm
Power input 200 VA
Power supply 230 V / 50 (60) Hz

Kneipp and **Hauff** bath

Electra® Hauff / Kneipp



CASCADE PLUS



Useful capacity 270 l
Water Pump 235 (375) l/min / 0,48 (0,17) MPa
Dimensions 1520 x 870 x 1035 mm
Power input 2400 (1500) VA
Power supply 3 x 400 V (230 V) / 50 (60) Hz

Whirling bath for **lower extremities**

Cascade®/ Plus



Useful capacity 30 l
Dimensions 950 x 870 x 1050 mm
Power input 800 VA
Power supply 230 V / 50(60) Hz

CORAL LYMFO



Useful capacity 30 l
Dimensions 950 x 870 x 1050 mm
Power input 850 VA
Power supply 3 x 400 V / 50 (60) Hz

Whirling bath for **upper extremities**

Coral®, Coral® Lymfo



Useful capacity 30 l
Dimensions 1060 x 660 x 795 - 1445 mm
Power input 600 VA
Power supply 230 V / 50 (60) Hz

LASTURA HOBBY



Useful capacity 40 l
Dimensions 1060 x 660 x 795 - 1445 mm
Power input 800 VA
Power supply 230 V / 50 (60) Hz

Whirling bath for **children / extremities**

Lastura® / Hobby



Useful capacity 54 l
Dimensions 880 x 760 x 790 mm
Power input 800 VA
Power supply 230 V / 50 (60) Hz

Whirling bath for **lower extremities**

Lastura®



Air consumption at 50 kPa 40 l / min
Weight 4 kg
No. of holes cca 2000

Bubble grid

Accessories